



Effective moisture management.




WICKS AWAY MOISTURE.

TruVal® is the only underpad on the market made with Capillary Surface Materials, a special class of novel cross-section yarns that dramatically enhance wicking performance. The extraordinary durability of the all-synthetic fabrics is enhanced even further by the laminated unipad construction, which also extends service life and reduces mechanical damage.

It is equipped with a napped barrier that prevents slipping and sliding, particularly on knit sheets. TruVal® folds more easily and efficiently, and lays flatter and more comfortably beneath the patient. It is available in two sizes, 33" x 36" and 36" x 41".

Kaiser Horizontal Wicking Rate Test



	TruVal®	lbex	Brushed Polyester
ZONE 1	4	4	4
ZONE 2	12	16	13
ZONE 3	34	186	123
ZONE 4	190	1055	2100

Values shown indicate time, in seconds, for 250 ml of liquid to reach outer edge of each zone.

For more information about Truval® incontinence products, or to contact a dedicated healthcare consultant, visit standardtextile.com or call 1.800.999.0400.