


# BARRIER SUPREME® REUSABLE SURGICAL WRAPPER



**BARRIER SUPREME®**  
Reusable Surgical Wrapper by Standard Textile

**STANDARD**  
TEXTILE est. 1940

# Reusable, cost-effective solution for wrapping surgical instrument sets, instrument trays and textile packs.



## FEATURES:

Barrier Supreme® Wrappers consist of a blended cotton/polyester to increase durability and resistance to punctures, runs, and snags. In addition, the blended fabric helps reduce the residual tape residue during the laundry process.

Standard Textile's Barrier Supreme® Surgical Wrappers have a liquid resistant finish. They feature color-coded edges to help simplify size identification and a QCM grid for monitoring laundering cycles.

## BEST PRACTICES:

To maximize performance, cost effectiveness, and longevity of reusable wrappers it is recommended to:

- Follow all appropriate AAMI and regulatory guidelines and instructions included in Standard Textile's IFU for Barrier Supreme® wrappers.
- When wrapping instrument sets and other metal items, always utilize corner guards to protect wrappers from damage due to sharp edges and corners.
- When securing the wrapper with tape, create a tab at the end of the tape by folding the end over on itself. This will make it easier to remove tape from the wrappers before sending the used wrapper to the laundry for washing and processing.
- If/When a wrapper is damaged and needs to be patched or repaired, be sure to follow AAMI guidelines for the inspection and patching/repair of reusable surgical textiles.

For more information about surgical wrappers, or to contact a dedicated healthcare consultant, visit [standardtextile.com](http://standardtextile.com) or call 1.800.999.0400.

**SurgiTex**<sup>™</sup>  
SURGICAL TEXTILE SOLUTIONS